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³ SIR (Society of Interventional Radiology). Guidelines for the use of retrievable and convertible vena cave filters: report from the Society of Interventional Radiology multidisciplinary consensus conference. J.A. Kaufman et al, Journal of Vascular Interventional Radiology 2006; 17 (3): 449-59.

⁴ Guidelines CIRSE 2007 (Cardiovascular and Interventional Radiological Society of Europe). Quality Improvement guidelines for percutaneous inferior vena cave filter placement for the prevention of pulmonary embolism.

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⁵ A prospective long-term study of 220 patients with a retrievable vena cava filter for secondary prevention of venous thromboembolism. P. Mismetti et al, en el pecho, 131 / 1, 2007.

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Tél : +33 (0)4 94 01 05 01

Site internet : <http://www.aln2b.com>

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INSTRUCTIONS FOR USE
VENA CAVA FILTER REMOVAL OR
REPOSITIONING KIT
JUGULAR APPROACH
(EN)
ALN 2 in 1 RS
Canada

First CE marking in 2010

Product reference: RK-2010

Ed. 03/2026

Contents

1: General information	33
2: Symbols	34
3: Kit content	35
4: ALN Vena Cava Filter extraction indications	36
5: Warnings.....	36
6: Potential Complications	41
7: ALN Vena Cava Filter extraction procedure.....	42
8: ALN Vena Cava Filter repositioning procedure	51
9: Patient's follow-up	55
10: Packaging and handling	56
11: Disposal of the device.....	56
12: Warranty.....	56

1: General information

The ALN Vena Cava Filter Removal or repositioning Kit is intended **only** to remove the ALN Vena Cava Filter with or without a hook when the indication of filtration no longer exists for a patient (see section 4. *ALN Vena Cava Filter extraction indications* and 5. *Warnings*), or to reposition the ALN Vena Cava Filter.

The clinical benefit of the ALN Vena Cava Filter Removal or repositioning Kit is to allow the removal of the ALN Vena Cava Filter when the indication for filtration no longer exists in a patient.

Its intended users are interventional radiologists, interventional cardiologists and vascular surgeons.

This medical device is available under three references:

- ALN Straight RS: reference FT.902010,
- ALN Pre-curved RS: reference FT.902010/VS2,
- ALN 2in1 RS: reference RK-2010.

The choice of device to use is individual to each patient. Nevertheless, ALN strongly recommends ensuring that both the Straight and Pre-curved references or the 2 in 1 reference are available.

The ALN Vena Cava Filter is an optional filter (definitive and / or temporary). Its efficacy and safety having been confirmed⁸, the ALN Vena Cava Filter can also be implanted permanently and there is no obligation to remove it.

The decision to remove the ALN Vena Cava Filter should be taken individually for each patient according to his/her clinical condition (benefit/risk ratio for the patient). Before making the decision to remove the filter, the degree of removal's urgency must be determined.

ALN 2 in 1 RS can only be used via the **jugular approach**.

Read the instructions carefully before use.

ALN 2 in 1 RS is supplied sterile (*ethylene oxide sterilization*). It is intended for one-time use and its components must not be re-sterilized.

⁸ A prospective long-term study of 220 patients with a retrievable vena cava filter for secondary prevention of venous thromboembolism. P. Mismetti et al, en el pecho, 131 / 1, 2007.

2: Symbols

This section presents the definitions of the symbols used on the product and its package.

	Catalog reference
	Batch number
	Serial number
	Use-by date
	Date of manufacture
	Warning
	Read instructions before using
	Sterilized using ethylene oxide, double sterile barrier system
	Do not reuse
	Do not resterilize
	Do not use if package is damaged
	Keep dry
	Keep away from sunlight
	Non-pyrogenic
	Medical Device
	Unique Device Identifier
	Manufacturer


3: Kit content

Ensure that all kit components are present before starting the procedure. If an item is missing, use another kit.

ALN 2 in 1 RS comprises the following items:

- A clamp catheter including:
 - a 7F sheath with two black marks (1),
 - a basket-shaped clamp comprised of eight 316LVM stainless steel strands (2) connected to a curved section (3). This curved section may be straight or curved, depending on the position of the introducer sheath. This curved section is connected by a central 316LVM stainless steel bar to a yellow cast-molded handle (4),
 - a three-way rinse valve (5),
- A dilator with red luer lock end piece (6),
- A 9F introducer sheath with radio-opaque marker (7),
- A puncture needle (8)
- A Teflon-coated J-guidewire. (9)

ALN 2 in 1 RS is provided sterile (*ethylene oxide sterilization*). It is intended for one-time use and its components must not be re-sterilized.

4: ALN Vena Cava Filter extraction indications

The indications for removal of the ALN Vena Cava Filter using ALN 2 in 1 RS listed below are approved by the SFICV⁹. They are based on the SIR¹⁰ and CIRSE¹¹ agreements:

- The permanent filtration indication is no longer present.
- The risk of symptomatic pulmonary embolism is acceptably low after successful initial anticoagulant treatment (therapeutic or preventive) or on change in the patient's clinical status.
- The patient is not likely to return to a high risk of pulmonary embolism following interruption of initial anticoagulant therapy and new clinical conditions.

In particular, do not remove the ALN Vena Cava Filter if the patient is to eventually undergo a high thromboembolic-risk surgical procedure.

- The patient's life expectancy is sufficiently long to achieve the estimated benefit from terminating filtration.

In particular, patients with a life expectancy of less than 6 months are not indicated for Vena Cava Filter removal.

- Extraction may be performed in a risk-free manner if the patient has no prior history of severe reaction to contrast agents, difficult venous approach, and if the extraction does not present significant risks for the inferior vena cava due to the nature, condition or position of the filter.
- The patient is aware of the potential risks associated with removal and agrees to the procedure being performed.
- The patient is aware that, in certain circumstances, the ALN Vena Cava Filter cannot be removed (technical or therapeutic impossibility) and accepts this risk.
- The patient is aware that the ALN Vena Cava Filter is an optional filter (permanent and/or temporary) and that it has been designed to be implanted permanently should the removal procedure be impossible.

5: Warnings

General warnings relative to the use of ALN 2 in 1 RS

1. Read the instructions for use carefully.
2. ALN 2 in 1 RS is intended **only**:
 - ✓ To remove the ALN Vena Cava Filter when the filtration indication no longer exists for a patient (see section 4. ALN Vena Cava Filter extraction indications).
 - ✓ Or to reposition the ALN Vena Cava Filter.

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¹⁰ SIR (Society of Interventional Radiology). Guidelines for the use of retrievable and convertible vena cave filters: report from the Society of Interventional Radiology multidisciplinary consensus conference. J.A. Kaufman et al, Journal of Vascular Interventional Radiology 2006; 17 (3): 449-59.

¹¹ Guidelines CIRSE 2007 (Cardiovascular and Interventional Radiological Society of Europe). Quality Improvement guidelines for percutaneous inferior vena cave filter placement for the prevention of pulmonary embolism.

3. ALN declines all liability in case of attempt to remove any other medical device using the ALN Vena Cava Filter Removal or repositioning Kit.
4. ALN declines all liability in case of attempt to remove the ALN Vena Cava Filter using a device other than the ALN Vena Cava Filter Removal or repositioning Kit provided for this purpose.
5. Any attempt to manipulate the ALN Vena Cava Filter using accessories other than those provided is prohibited and under the practitioner's entire responsibility. Except a snare with fitting dimensions that can be used to remove an ALN Vena Cava Filter with hook.
6. ALN 2 in 1 RS is provided ready to use, sterile (*ethylene oxide sterilized*) and non-pyrogenic. Do not use the product if its package is damaged or open, or if there is a doubt whether is sterile.
7. The product should be stored in a dry place and kept away from sunlight.
8. Do not use ALN 2 in 1 RS beyond the expiration date indicated on the label.
9. ALN 2 in 1 RS is intended for one-time use. It must not be reused and the kit's components must not be re-sterilized. Reuse even after a new cleaning and a new sterilization could expose the patient to an infectious risk and a potential malfunction of the device (device not controlled by the manufacturer).
10. ALN 2 in 1 RS can **only** be used via the **jugular approach**.
11. The ALN Vena Cava Filter is an optional filter (definitive and / or temporary). Its efficacy and safety having been confirmed¹², the ALN Vena Cava Filter can also be implanted permanently and there is no obligation to remove it.
12. The decision to remove the ALN Vena Cava Filter should be taken individually for each patient according to his / her clinical condition (benefit / risk ratio for the patient). Before making the decision to remove the filter, the degree of removal's urgency must be determined.

General warnings relative to the patient's clinical condition

13. Patients with still active thromboembolic disease undergoing anticoagulant treatment:
 - ✓ The time during which the patient is exposed to the risk of relapse of pulmonary embolism must be individualized to each patient (integration of known risk factors).
 - ✓ The anticoagulation of a patient for whom removal of the ALN Vena Cava Filter is planned must be taken into account or monitored/evaluated before the intervention so as to avoid any hemorrhagic risk.
 - ✓ Any suspicion of failure or complication of the anticoagulant treatment must be eliminated prior to considering removal of the ALN Vena Cava Filter.
14. Patients with no thromboembolic disease:
 - ✓ Removal of the ALN Vena Cava Filter, for patients at risk of thromboembolic disease, with no diagnosis of thromboembolic disease, for whom a filter has been implanted for preventive purposes, may be considered once the initial prevention has been successfully completed and is likely to be maintained, or once the high risk of pulmonary embolism has been reduced.
 - ✓ There should be no obvious evidence of development of thromboembolic disease before removing the ALN Vena Cava Filter.

¹² A prospective long-term study of 220 patients with a retrievable vena cava filter for secondary prevention of venous thromboembolism. P. Mismetti et al, en el pecho, 131 / 1, 2007.

15. Other cases:

- ✓ If, for specific reasons, the ALN Vena Cava Filter is no longer able to play its preventive role with respect to pulmonary embolism, or if it represents a hazard to the patient, removal must be considered.
- ✓ If the filtration indication is still present, another ALN Vena Cava Filter should be implanted.

Warning before considering the removal of the ALN Vena Cava Filter

16. Control imaging and a precise blood evaluation must be performed.

In all cases, before any ALN Vena Cava Filter removal procedure, check the permeability of the inferior vena cava and ensure there are no thrombi in the vicinity of the ALN Vena Cava Filter. This examination can be performed via Doppler ultrasound, CT scan with injection, or cavography.

If a thrombus is found to be trapped in the ALN Vena Cava Filter, its size and type must be determined prior to making the removal decision. Indeed, the presence of a thrombus trapped in the ALN Vena Cava Filter may indicate that the filter has played its role (in this case, ensure that the disease is not progressive), or that the ALN Vena Cava Filter presents a thrombosis problem.

Continuation or termination of anticoagulant treatment should be determined on a per case basis by the medical team.

Cavography-related warnings

17. Anterior posterior (AP) and lateral cavography must be performed before the removal procedure. The purpose of this examination is to check:

- ✓ the permeability and anatomy of the inferior vena cava,
- ✓ the absence of thrombi,
- ✓ the position of the ALN Vena Cava Filter.

The result of the cavography in terms of inferior vena cava anatomy, along with the position of the ALN Vena Cava Filter, shall assist the practitioner in selecting the type of Vena Cava Filter Removal or repositioning Kit to use (Straight - reference FT.902010, Pre-curved - reference FT.902010/VS2, or 2 in 1 - reference RK-2010).

18. During injection of contrast medium into the introducer sheath, do not exceed a pressure of 450 psi.
19. Ensure that the cavography can be performed without risk for the patient: take into account any problems associated with use of contrast media (allergy, renal failure).
20. If there is any doubt concerning the anatomy of the vena cava and / or the position of the ALN Vena Cava Filter, new cavography should be taken with different incident angles.
21. A control cavography may be performed after the extraction procedure.
22. If the ALN introduction system (introducer sheath and dilator) is used to perform a control cavography, **ensure that the dilator has been removed before performing the cavography**. Not removing the dilator could damage the vena cava.
23. If the control cavography is performed using the introducer sheath of ALN 2 in 1 RS, ensure that the distal end of this sheath does not enter into contact with the implanted ALN Vena Cava Filter. This could modify the filter's position (risk of tilting or migration).

General procedure-related warnings

24. The expertise of physician and non-physician operators who will be performing the procedures must conform to the recommendations of the **SFICV**¹³.
- ✓ A junior operator must have previously assisted a senior during at least 6 removals (as 2nd operator) and must have removed at least 6 filters himself/herself under the supervision of a senior to be deemed competent.
 - ✓ A senior may maintain his/her skill by removing at least 6 filters per year without incident.
 - ✓ The electro radiology operator should possess basic vascular procedure training, as defined by the **SFICV**⁶.
 - ✓ Nurses, including operating theatre nurses, must have undergone specific invasive and interventional vascular imaging training.
25. ALN warns users of the increased risk of migration in the following cases:
- ✓ major surgery (thorax, abdomen, etc.),
 - ✓ procedures potentially modifying haemodynamics,
 - ✓ procedures that may cause contact with the inferior vena cava and/or the filter implantation site.

Being aware of the above-mentioned potential risks, it is important to ensure that all elements pertaining to ALN Vena Cava Filter implantation are reported in the patient record, which should be checked prior to each procedure.

26. Operators must wear sterile clothing.
27. Before the procedure, the elements of the kit (J-guidewire and introduction system (dilator + introducer sheath)) should be rinsed with a sterile physiological saline solution or similar isotonic solution.
28. The entire procedure must be performed under medical imaging (preferably fluoroscopy).

Warnings relative to use of the J-guidewire

29. The J-guidewire should be used by a physician with appropriate cardiac catheterization training. Never advance the J-guidewire or introduction system without radiological control.
30. During J-guidewire progression, check for the absence of bends or loops under medical imaging. Always ensure that minimum 2-3cm of the J-guidewire extend beyond the proximal end of the introduction system.
31. Do not attempt to advance or remove the J-guidewire in case of resistance. The cause of resistance should be determined by radioscopy.
32. Make sure to remove the J-guidewire and the dilator before performing the cavography in order to avoid any trauma to the vena cava.
33. Radiographies performed without radio contrast medium, that do not clearly show the anatomy of the inferior vena cava, may be misinterpreted.

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Warnings relative to implementation of the introducer sheath

34. Ensure that the distal end of the introducer sheath does not enter into contact with the implanted ALN Vena Cava Filter. This could modify the filter's position (risk of tilting or migration).
35. Before inserting the clamp catheter into the introducer sheath, the radio-opaque band on the sheath's distal end must be positioned 5 cm above the filter's apex. This distance is essential to start the procedure under optimal conditions.

Warnings relative to insertion of the clamp catheter into the introducer sheath.

36. Ensure that the J-guidewire and the dilator have been removed before inserting the clamp catheter.
37. To avoid kinking the introducer sheath, ensure that the introducer sheath is handled by its base and is not bent.
38. Always ensure that the introducer sheath is adequately rinsed with a sterile heparinized physiological saline solution or similar isotonic solution. Insufficient rinsing of the introducer sheath during the procedure may cause clot formation in it.
39. Never insert the clamp catheter into the introducer sheath before having retracted the clamp of the clamp catheter into its catheter. For this, move back the clamp inside its catheter by rotating the yellow luer lock left to unlock, then pulling it for approximately 5 cm.
40. Ensure that the base of the introducer sheath is in contact with the base of the clamp catheter.
41. Use the yellow handle luer lock and never unlock the luer lock screwed onto the clamp catheter sheath.
42. Lock the clamp catheter luer lock by rotating the yellow handle right, causing the clamp to open.

Warnings relative to capturing the ALN Vena Cava Filter

43. Never attempt to remove the ALN Vena Cava Filter by one of its legs. The ALN Vena Cava Filter must strictly be removed by its apex.
44. To remove the ALN Vena Cava Filter under optimal conditions, the filter's apex must be fully engaged into the clamp of the clamp catheter.
45. Depending on the position of the ALN Vena Cava Filter inside the vena cava, it may be preferable to grasp the filter's apex with the curved section in straight or curved position. The angle of the curved section can be modified by lowering the introducer sheath on the curved section. When the base of the introducer sheath is fully in contact with the clamp catheter base, the clamp of ALN 2 in 1 RS is curved in its maximum. When the introducer sheath is lowered by about 1,5 cm, the clamp of ALN 2 in 1 RS is getting straight.
46. When the filter is ready to be captured, do not pull on the clamp catheter until the filter is completely captured.
47. If the position of the filter does not allow easy removal (filter' apex against the wall of the vena cava, for example), the degree of urgency and the operator's experience should be taken into account before performing any further endovascular maneuvers.
48. Losing the filter during the procedure being a possible risk, the recommendations made by ALN must be followed.

Filter recovery-related warnings

49. Do not pull on the clamp catheter until the first black mark on the clamp catheter is not visible or has been passed.

50. Check the progression of the filter and the clamp catheter in the introducer sheath under medical imaging
51. In order to remove the ALN Vena Cava Filter from the introducer sheath and to use the sheath for postoperative cavography (recommended), sufficient tension must be applied to the filter.
52. Incorrect handling of the ALN vena cava filter may expose the patient to a risk of sting/injury or infection.

Warning after the removal of the ALN Vena Cava Filter

53. A medical imaging is essential after the removal of the ALN Vena Cava Filter in order to confirm the absence of inferior vena cava injuries, or the presence of a thrombus in the vena cava.

If any inferior vena cava injury is observed, some additional procedures and close clinical monitoring will be necessary.

54. The integrity of the ALN Vena Cava Filter must be confirmed.

Should a part of the filter be missing, check the sheath and check the abdomen with radiology images to locate and document the missing filter elements.

55. Therapeutic follow-up after the ALN Vena Cava Filter removal.

Patients should be treated according to the diagnosis of their thromboembolic disease, with no specific treatment and continuing their therapy if applicable.

All patients should be followed up for a new, recurrent, or progressive deep vein thrombosis and/or pulmonary embolism.

Other warnings

56. After use, the kit components should be discarded according to usual medical practice and applicable regulations.

6: Potential Complications

Procedural complications reported with filter retrieval include (but are not limited to):

- Migration of filter fractured components
- Thrombus migration during filter retrieval
- Cardiac tamponade after filter removal
- Caval injury such as laceration, development of a pseudoaneurysm, rupture and intussusception,
- Pulmonary embolism,
- Deep vein thrombosis,
- Caval collapse at the point of prior prolonged filter endothelialization,
- Acute caval thrombus.

A certain number of potential J-guidewire-related complications exist when used on patients with Vena Cava Filters. These complications include, in particular, a risk of migration, of tilting and the possibility of the J-guidewire becoming trapped in the filter¹⁴. Consequently, the extraction procedure must be performed by practitioners with appropriate cardiac catheterization training.

¹⁴ Guide Wire Entrapment by Inferior Vena Cava filters: an experimental study. M.J. Rosen et al, Journal of American College of Surgeons, Vol 201, No3, 2005.

7: ALN Vena Cava Filter extraction procedure

Warning: The entire procedure must be conducted under medical imaging.

Warning: Before starting the ALN Vena Cava Filter removal procedure, several checks need to be performed. See the section "**Warning before considering the removal of the ALN Vena Cava Filter**" in this manual. In particular, a medical imaging is essential in order to check the presence of a thrombus in the ALN Vena Cava Filter. If a thrombus is present, do not extract the ALN Vena Cava Filter.

Warning: ALN 2 in 1 RS can only be used via the **jugular approach**.

1: Patient preparation:

Disinfect and position sterile tissues (1). The necessary equipment is carefully prepared (2) (extraction kit, sterile heparinised physiological saline solution or similar isotonic solution for rinsing the system, etc.).

(1)

(2)

After local anesthetic, puncture the right internal jugular vein using the provided puncture needle (3).



(3)

2: Inserting the J-guidewire and rinsing the introduction system

Warning: The J-guidewire should be used by a physician with appropriate cardiac catheterization training. Never advance the J-guidewire or the introduction system without radiological control.


Flush the J-guidewire with the sterile physiological saline or similar isotonic solution, fixing a syringe to the luer lock base of the J-guidewire's distributor (4). Insert the flexible end of the J-guidewire (5) into the puncture needle.



Under medical imaging, introduce the J-guidewire through the needle using the stiffener and advance the J-guidewire until it is positioned above the filter.

Hold the J-guidewire in place and remove the puncture needle and the J-guidewire's stiffener (6).

To facilitate the progression of the introduction system, we recommend making a sufficiently large incision to the tissues surrounding the J-guidewire.

***Warning:** During the J-guidewire progression, check for the absence of bends or loops under medical imaging. Always ensure that 2-3 cm of the J-guidewire extends beyond the proximal end of the introduction system (7).*


***Warning:** Do not attempt to advance or remove the J-guidewire in case of resistance. The cause of resistance should be determined by radioscopy. If the cause of the resistance cannot be determined, remove the J-guidewire and the introduction system.*

3: Advancing the introduction system (introducer sheath + dilator)


Under medical imaging, slide the introduction system (8) over the J-guidewire and position the radio-opaque band 5 cm above the filter's apex (9). This distance is essential to start the procedure under optimal conditions. If this position is not respected, you may open the extraction clamp inside the filter (10).

(8)

(9)


(10)

Warning: Ensure that the distal end of the dilator does not enter into contact with the implanted ALN Vena Cava Filter. This could modify the ALN Vena Cava Filter's position (risk of tilting or migration).


Warning: To avoid kinking the introducer sheath, ensure that the introducer sheath is handled by its base and is not bent.

Remove the J-guidewire. (11)

(11)

Remove the dilator by rotating the red luer lock counter-clockwise. (12)

(12)

4: Inserting the clamp catheter into the introducer sheath

Warning: Ensure that the J-guidewire and the dilator have been removed before inserting the clamp catheter.

Before inserting the clamp catheter into the introducer sheath, rinse the sheath and rinse the clamp catheter using the 3-way valve with heparinized physiological saline.

Warning: Insufficient rinsing of the introducer sheath during the procedure may cause clot formation in it. Move back the clamp inside its 7F sheath by rotating the yellow luer lock left to unlock, then pulling it for approximately 5 cm (13).

(13)

Warning: Never insert the clamp catheter into the introducer sheath before having retracted the of the clamp catheter into its catheter.

Warning: Use the yellow handle luer lock and never unlock the luer lock screwed onto the clamp catheter (14).



Luer lock of the yellow handle

Luer lock screwed onto the clamp catheter

Under frontal and lateral imaging, fully insert the clamp catheter into the introducer sheath (15).

(15)

Warning: Ensure that the base of the introducer sheath is in contact with the base of the clamp catheter. (16)

(16)

Push, then lock the clamp catheter's yellow luer lock handle, rotating it to the right (17) in order to release the curved section, which is then in curved position, along with the clamp's hooks.

(17)

Depending on the position of the ALN Vena Cava Filter inside the vena cava, it may be preferable to grasp the filter's apex with the curved section in straight or curved position. The angle of the curved section can be modified by lowering the introducer sheath on the curved section. When the base of the introducer sheath is fully in contact with the clamp catheter base, the clamp of ALN 2 in 1 RS is curved in its maximum (18).

(18)

When the introducer sheath is lowered by about 1 cm, the angle of the curved section decreases. (19)

(19)

When the introducer sheath is lowered by about 1,5 cm, ALN 2 in 1 RS is getting straight (20).

(20)

The clamp can be rotated by turning ALN 2 in 1 RS handle.

5: Capturing the ALN Vena Cava Filter

Warning: Never attempt to remove the ALN Vena Cava Filter by one of its legs. The ALN Vena Cava Filter must strictly be removed by its apex. (21)

(21)

Slowly move down the clamp catheter over the apex of the filter and ensure that the clamp catheter is fully in contact with the apex of the ALN Vena Cava Filter (22) (under medical imaging, using preferably a lateral or a three-quarter angle).

(22)

Warning: Depending on the position of the ALN Vena Cava Filter inside the vena cava, it may be preferable to grasp the filter's apex with the curved section in straight or curved position. The angle of the curved section can be modified by lowering the introducer sheath on the curved section. When the base of the introducer sheath is fully in contact with the clamp catheter base, the clamp of ALN 2 in 1 RS is curved in its maximum. When the introducer sheath is lowered by about 1,5 cm, the clamp of ALN 2 in 1 RS is getting straight.

When the base of the introducer sheath is fully in contact with the clamp catheter base, the clamp is curved in its maximum. If the ALN Vena Cava Filter is centered inside the vena cava, lower the introducer sheath by 1,5 cm to straighten the curved section. ALN 2 in 1 RS is now straight. To grasp the filter's apex, slowly lower the clamp until it makes contact with the apex of the ALN Vena Cava Filter (23).

(23)

When the base of the introducer sheath is fully in contact with the clamp catheter base, the clamp is curved in its maximum. If the ALN Vena Cava Filter is tilted inside the vena cava, orient the clamp using the handle such that ALN 2 in 1 RS is pointing at the apex of the ALN Vena Cava Filter. Next, slowly lower the clamp onto the apex of the ALN Vena Cava Filter until it is in contact with the apex. (24)

(24)

Warning: If the position of the ALN Vena Cava Filter does not allow an easy removal (filter's apex against the wall of the vena cava for example), the degree of urgency and the operator's experience should be taken into account before performing any further endovascular maneuvers.


First, move down the introducer sheath slowly over the ALN Vena Cava Filter until the black mark on the clamp catheter is visible. This mark indicates that the filter has entered into the introducer sheath (25).

(25)

Warning: Do not pull on the clamp catheter until the first black mark on the clamp catheter is not visible or has been passed.

If you want to use the introducer sheath to perform a control cavography after the ALN Vena Cava Filter removal, pull back the clamp catheter with a slow and steady movement until the second black mark appears (26). Then, pull back firmly to remove the clamp catheter and the filter from the introducer sheath (27).

(26)


(27)

Warning: Check the progression of the filter and the clamp catheter in the introducer sheath under medical imaging

Take out the ALN Vena Cava Filter from the introducer sheath. (28)

Warning: In order to remove the ALN Vena Cava Filter from the introducer sheath and to use the sheath for postoperative cavography (recommended), sufficient tension must be applied to the ALN Vena Cava Filter.

(28)

Perform a cavography using the introducer sheath left in place.

6: End of procedure

Remove the introducer sheath from the inferior vena cava and achieve hemostasis at the puncture site.

8: ALN Vena Cava Filter repositioning procedure

We recommend repositioning the ALN Vena Cava Filter in the following cases:

- Premature ALN Vena Cava Filter deployment during placement,
- ALN Vena Cava Filter deployment in a position different to that initially planned,
- Significant tilting reducing the ALN Vena Cava Filter's filtration capacity.

1: Removing the ALN Vena Cava Filter

Slowly move down the clamp catheter over the filter's apex and ensure that the clamp catheter is fully in contact with the apex of the ALN Vena Cava Filter (29) (under medical imaging, using preferably a lateral or a three-quarter angle).

(29)

Warning: Depending on the position of the ALN Vena Cava Filter inside the vena cava, it may be preferable to grasp the filter's apex with the curved section in straight or curved position. The angle of the curved section can be modified by lowering the introducer sheath on the curved section. When the base of the introducer sheath is fully in contact with the clamp catheter base, the clamp of ALN 2 in 1 RS is curved in its maximum. When the introducer sheath is lowered by about 1,5 cm, the clamp of ALN 2 in 1 RS is getting straight.

When the base of the introducer sheath is fully in contact with the clamp catheter base, the clamp is curved in its maximum. If the ALN Vena Cava Filter is centered inside the vena cava, lower the introducer sheath by 1,5 cm to straighten the curved section. ALN 2 in 1 RS is now straight. To grasp the filter's apex, slowly lower the clamp until it makes contact with the apex of the ALN Vena Cava Filter. (30)

(30)

When the base of the introducer sheath is fully in contact with the clamp catheter base, the clamp is curved in its maximum. If the ALN Vena Cava Filter is tilted inside the vena cava, orient the clamp using the handle such that the clamp of ALN 2 in 1 RS is pointing at the apex of the ALN Vena Cava Filter. Next, slowly lower the clamp onto the apex of the ALN Vena Cava Filter until it is in contact with the apex. (31)

(31)

Warning: If the position of the ALN Vena Cava Filter does not allow an easy removal (filter's apex against the wall of the vena cava, for example), the degree of urgency and the operator's experience should be taken into account before performing any further endovascular maneuvers.

In order to reposition the ALN Vena Cava Filter, slowly move down the introducer sheath over the ALN Vena Cava Filter until the black mark on the clamp catheter is visible. This mark indicates that the filter has entered into the introducer sheath. (32)


(32)

Warning: At this step of the procedure, never pull on the clamp catheter.


Warning: Once the ALN Vena Cava Filter is inside the introducer sheath, a solution of sterile physiological saline must be injected regularly in order to avoid the formation of blood clots that could not allow deploying the filter correctly.

2: Positioning the introducer sheath at the implantation site

Position the radio-opaque band of the introducer sheath approximately 55 mm for the ALN Vena Cava Filter without hook (33.1) or 60 mm for the ALN Vena Cava Filter with hook (33.2) below the lowest renal vein if the possibility to implant a filter under the lowest renal vein has been confirmed.

**(33.1) ALN Vena Cava Filter
without hook**


**(33.2) ALN Vena Cava Filter
with hook**

Warning: Do not release the filter as long as its position is not optimal, the most appropriate implantation site should be confirmed by cavography.

Warning: Suprarenal implantation of the ALN Vena Cava Filter should be only performed after checking the diameter's size of the suprarenal vena cava (maximum 32 mm) and the room available between the highest renal vein and the atrium, at least 55 mm (height of the filter without hook) or 60 mm ((height of the filter with hook).



3: Filter release at the implantation site

Holding the clamp handle, withdraw the introducer sheath back over the clamp catheter with a slow and steadily retrieval movement until the clamp is free. The clamp should now be fully open, and the filter completely deployed. (34)


(34)

Warning: Do not pull up too much the introducer sheath in order to avoid opening the curved section. Actually, if the curved section is open, ALN 2 in 1 RS will be with an angulation and the filter could not be implanted in an optimal way in the vena cava. (35)

(35)

Next pull back the handle with a slow and straight movement in order to get the clamp away from the filter's apex. The clamp should be 2 cm away from the filter's apex. (36) The introducer sheath and the handle should be pulled back at the same time.

(36)


Warning: The retrieval movement should be slow and straight to avoid any tilting of the filter.

Warning: Pull the clamp catheter with the introducer sheath, otherwise the clamp may be closed over the apex of the filter.

Warning: A misplaced filter which nevertheless protects the patient against a pulmonary embolism in an appropriate way should be left in place. If a misplaced filter, doesn't protect the patient against a pulmonary embolism in an appropriate way, it should be repositioned, or a second filter should be implanted. Removal, either by surgical or by using one of the ALN Vena Cava Filter Removal or repositioning Kit intended for this purpose, is recommended for a misplaced filter.

4: Removing the clamp catheter

In order to remove the clamp catheter, slowly move down the introducer sheath onto the clamp catheter until it is completely covered. (37)

(37)

Warning: The procedure must be followed under medical imaging in order to ensure that the clamp is fully inserted into the introducer sheath.

Withdraw the clamp catheter by a slowly retrieval movement until the clamp goes out from the base of the introducer sheath.

Warning: Check the progression of the clamp catheter under medical imaging.

Perform a cavography using the introducer sheath left in place.

5: End of procedure

Remove the introducer sheath from the inferior vena cava and achieve hemostasis at the puncture site.

9: Patient's follow-up

Considering the above-mentioned risk, patient's follow-up should include the following points:

- ✓ Examination of the catheter insertion site to make sure about a good healing.
- ✓ A medical imaging is essential after the removal of the ALN Vena Cava Filter in order to confirm the absence of inferior vena cava injuries, or the presence of a thrombus in the vena cava.

This imaging is even more essential if the procedure was long and difficult, or if the patient complained about pain during or after this procedure.

If any inferior vena cava injury is observed, some additional procedures and close clinical monitoring will be necessary.

- ✓ The integrity of the ALN Vena Cava Filter must be confirmed.

Should a part of the filter be missing, check the sheath and check the abdomen with radiology images to locate and document the missing filter elements.

An additional procedure is possible by the decision of the physician who conducted the initial procedure.

- ✓ Therapeutic follow-up after the removal of the ALN Vena Cava Filter.

Patients should be treated according to the diagnosis of their thromboembolic disease, with no specific treatment and continuing their therapy if applicable.

All patients should be followed up for of a new, recurrent, or progressive deep vein thrombosis and/or pulmonary embolism.

The resumption of an anticoagulant treatment must be adapted by the medical team to each patient.

A Summary of Safety and Performance Characteristics will be available in the European Medical Device Database (EUDAMED) via this link: [SS\(C\)P > SS\(C\)P: DOC-226 - EUDAMED](#).

The public website of Eudamed is: <https://ec.europa.eu/tools/eudamed/#/screen/home>.

An electronic copy of this user manual is available on the ALN website: www.aln2b.com.

10: Packaging and handling

ALN 2 in 1 RS is provided ready to use, sterile (ethylene oxide sterilized) and non-pyrogenic. Do not use the product if its packaging is damaged or open, or if there is a doubt whether is sterile.

The product should be stored in a dry place and kept away from sunlight.

Do not use ALN 2 in 1 RS beyond the expiration date indicated on the label.

ALN 2 in 1 RS is intended for one-time use. It must not be reused and the kits components must not be re-sterilized.

After use, the kit components should be discarded as per usual medical practice and according to applicable regulations.

11: Disposal of the device

A used device becomes infectious risk healthcare waste. It must undergo decontamination and be eliminated by the appropriate channel (example: incineration or pre-treatment by disinfection) in accordance with usual medical practices and according to the regulations in force.

An unused device must be disposed of as follows: the needle must be disposed of with waste from health care activities at risk of infection if it has been removed from its packaging because it may expose to a risk of injury. The cardboard box and the paper instructions are to be recycled. The rest of the device can be disposed of through the conventional waste stream of healthcare establishments.

12: Warranty

ALN declares to have manufactured this product with all necessary care. This guaranty replaces and excludes all other formal or tacit guaranties granted by virtue of any other provision not expressly stated hereto, in particular all tacit guaranties concerning the product's suitability for sale and for specific use, due to handling, storage, cleaning and sterilisation of this product, along with certain factors related to the patient, diagnosis, treatment and surgical procedures, etc., beyond the control of ALN, have an impact on this product and on the results achieved with it.

ALN cannot be held liable for any fortuitous or indirect damage, nor for losses or costs resulting, either directly or indirectly, from the use of this product. It will, however, be obliged to replace it.

Furthermore, ALN accepts no other additional commitments or liabilities in relation to this product and does not authorize any third parties to do so.

FDA Approval

ALN S.A.R.L.
589, Chemin du Niel
83230 Bormes les Mimosas
France

Tel: +33 (0)4 94 01 05 01

Website: <http://www.aln2b.com>
Email address: contact@aln2b.com